UK

PHYSICIAN
INSTRUCTIONS

DuraPlug®
Synthetic Extended Temporary
Punctal/Canalicular Insert

DESCRIPTION
DuraPlug Synthetic Extended Temporary Punctal/Canalicular Inserts are intended for temporary use with patients experiencing dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations which can be relieved by blockage of the canaliculus. The DuraPlug Synthetic Extended Temporary Punctal/Canalicular Insert is made from synthetic absorbable suture material, E-Caprolactone-L-Lactide copolymer (PCL). The plug is provided dyed violet with D&C Violet No. 2 and is coated. The substances contained in the coating are noncollagenous and nonantigenic.

The plugs are available in one length, 2.0mm, and in three diameters:

0.2mm, 0.3mm, and 0.4mm. There are two plugs per package so that a plug can be placed in each canaliculus.

INTENDED USE
DuraPlug Synthetic Extended Temporary Punctal/Canalicular Inserts are designed to be inserted through the punctal opening into the canaliculus to block tear drainage through the lacrimal drainage system.

INDICATIONS
DuraPlug Synthetic Extended Temporary Punctal/Canalicular Inserts may be used in the treatment of Dry Eye Syndrome and the dry eye components of various ocular surface diseases. When indicated, DuraPlug Synthetic Extended Temporary Punctal/Canalicular Inserts may be used after ocular surgery to prevent complications due to dry eye and to enhance the retention of ocular medications. Patients experiencing dry eye related contact lens problems also may be aided by DuraPlug Synthetic Extended Temporary Punctal/Canalicular Inserts.

CONTRAINDICATIONS
DuraPlug Synthetic Extended Temporary Punctal/Canalicular Inserts are contraindicated for patients who are experiencing epiphora, inflammation of the eyelid, or tearing secondary to dacryocystitis with mucopurulent discharge. If the patient experiences irritation or epiphora after insertion of DuraPlug Synthetic Extended Temporary Punctal/Canalicular Inserts, a saline irrigation or a probe may be used to expel the plug through the lacrimal sac into the nose or throat.

PATIENT PREPARATION
Inspect the patient's punctum to determine the appropriate size plug prior to insertion of the DuraPlug Synthetic Extended Temporary Punctal/Canalicular Insert.

DIRECTIONS FOR USE
Remove the appropriate diameter plug from the plug holder using jewelers' forceps. While holding the plug with the forceps, place the plug partially through the punctal opening. Release the plug and use the tip of

the forceps to push the plug the rest of the way through the punctal opening and into the canaliculus. A small amount of lubricant may be used on the plug in order to facilitate insertion.

ADDITIONAL INFORMATION
Effective occlusion in the first days after insertion is dependent on the size of the patient's punctal opening and the size plugs used. Absorption of DuraPlug Synthetic Extended Temporary Punctal/Canalicular Insert is essentially completed between 60 and 180 days. The patient may experience rapid improvement of dry eye symptoms during the initial days followed by a return to previous discomfort. Such signs of improvement suggest that the patient may benefit by permanent, yet reversible, punctal occlusion with EagleVision Silicone Punctal Plugs.

CAUTIONS
US Federal Law restricts this device for sale by or on the order of a medical physician. Care should be used not to perforate the punctum or

canaliculus while inserting the DuraPlug Synthetic Extended Temporary Punctal/Canalicular Insert. Perforation may cause pain and increase the risk of infection. If perforation occurs, delay the placement of the plug until the wound heals. DuraPlug Synthetic Extended Temporary Punctal/Canalicular Inserts may enhance the effect of ocular medication used on the eye. Depending on the type of medications being used, dosage may need to be reduced accordingly.

PRECAUTIONS
Contents are sterile if individual package is unopened and undamaged. Do not use if package has been previously opened.

STORAGE
Ambient conditions apply. Do not expose to extreme temperatures.

SYMBOLS USED ON LABELING

- 

* Do not reuse
- 

* Use until month & year

STERILEEO

* Sterile unless package is opened or damaged.
Method of Sterilization – Ethylene Oxide.

LOT

* Batch Number

* See Instructions for Use

ECREP

European Representative:
Quality First International
11 Burford Road, Suites 317-318
Stratford, London E15 2ST
UK
Telephone: 44 208 221 2361
Telefax: 44 208 221 1912

Manufactured for:
EagleVision, Inc.
8500 Wolf Lake Drive, Suite 110
Memphis, TN 38133 USA

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